

Product Datasheet

MAP LC3 Beta antibody (orb157847)

Description	MAP LC3 Beta antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Microtubule-associated proteins 1A/1B light chain 3A (25-121/121 aa)
Target	MAP1LC3A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14 kDa
Uniprot ID	Q9H492
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.