



Product Datasheet

MC5 Receptor antibody (orb157857)

Description	MC5 Receptor antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human MC5 Receptor (78-125/325 aa)
Target	MC5 Receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P33032
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse pl...



Flow cytometric analysis of spleen cell

...