


Product Datasheet

MGLUR3 antibody (orb157870)

Description	MGLUR3 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Metabotropic Glutamate Receptor 3 (365-460/879 aa)
Target	MGLUR3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	96 kDa
Uniprot ID	Q14832
Dilution Range	WB=1:500-2000, Flow-Cyt=5µg/Test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Mouse Liver

...