



Product Datasheet

MGST1 antibody (orb157871)


Description	MGST1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MGST1 (81-155/155 aa)
Target	MGST1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17 kDa
Uniprot ID	P10620
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat stom...



Immunohistochemical staining of human ki...



Immunohistochemical staining of human ki...