



Product Datasheet

MUC5AC antibody (orb157913)


Description	MUC5AC antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MUC5AC/Mucin 5AC (601-700/5030 aa)
Target	MUC5AC
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	523 kDa
Uniprot ID	P98088
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human
stomach
using MU...



IHC-P
analysis of
rat stomach
using MUC5...



IHC-P
analysis of
mouse
stomach
tissue u...