



Product Datasheet

NDUFS7 antibody (orb157939)

Description	NDUFS7 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human NDUFS7 (101-160/213 aa)
Target	NDUFS7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	075251
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse skeletal muscle tissue.



Western blot analysis of mouse muscle and heart tissue.