



Product Datasheet

NLK (Phospho-Thr298) antibody (orb157946)

Description	NLK (Phospho-Thr298) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Nemo-like kinase around the phosphorylation site of Thr298 (HM(p-T)QE)
Target	Phospho-NLK (Thr298)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	Q9UBE8
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Sample: Cerebral cortex (Mouse)
Lysate a...



Paraformaldehyde-fixed, paraffin embedde...