


Product Datasheet

CHPT1 antibody (orb1580)

Description	CHPT1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human CHPT1 (61-160/406 aa)
Target	CHPT1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	Q8WUD6
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Cerebrum (Mouse);Heart (Mou...