

Product Datasheet

phospho-GSK-3 Beta (Ser9) Antibody (Biotin) [Out of stock] (orb1580356)

Description

phospho-GSK-3 Beta (Ser9)

Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GSK-3 Beta around the phosphorylation site of Ser9
Target	phospho-GSK-3 Beta (Ser9)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P49841
Dilution Range	WB: 1:500-2000, IHC-P: 1:20-200, IHC-F: 1:20-200, ICC: 1:20-200, IF: 1:20-200
Expiration Date	12 months from date of receipt.