

Product Datasheet

Phospho-Smad2 (Ser465 + Ser467) Antibody (BF350) (orb1581123)

Description	Phospho-Smad2 (Ser465 + Ser467) Antibody (BF350)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF350
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SMAD2 around the phosphorylation site of Ser465/467
Target	Phospho-Smad2 (Ser465 + Ser467)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15796
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.