



Product Datasheet

PDL1 antibody (orb158130)


Description	PDL1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human CD274 (64-160/290 aa)
Target	PD-L1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	Q9NZQ7
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Flow cytometry of
Mouse spleen
cells(blu...



Flow cytometry of
Hela.Primary
Antibody ...



Western
blot:Du145(Human)
Cell;Lovo(Huma...