



Product Datasheet

PHEX antibody (orb158144)


Description	PHEX antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PHEX. Please contact us for the exact immunogen sequence. The peptide is available as orb374817.
Target	PHEX
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	-100 µg (in 200 µl): 0.5 mg/ml-200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86 kDa
Entrez	5251
Dilution Range	WB: 1:100-1000, IHC-P: 1:100-600 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat panc...



WB analysis of rat pancreas (lane 1); ra...



IHC-P image of mouse ovary tissue using ...