

Product Datasheet

PIK3 gamma (Phospho-Ser1101) antibody (orb158152)

Description	PIK3 gamma (Phospho-Ser1101) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PIK3 gamma around the phosphorylation site of Ser1101 (KH(p-S)A)
Target	phospho-PIK3 gamma (Ser1101)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	121 kDa
Uniprot ID	P48736
Expiration Date	12 months from date of receipt.