


Product Datasheet

PIK3 gamma (Phospho-Thr1024) antibody (orb158153)

Description	PIK3 gamma (Phospho-Thr1024) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PIK3 gamma around the phosphorylation site of Thr1024 (HH(p-T)NL)
Target	phospho-PIK3 gamma (Thr1024)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	121 kDa
Uniprot ID	P48736
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Jurkat (Human) Cell; NIH/3T3 (...)