

Product Datasheet

CHRM4 Antibody (BF350) (orb1582144)

Description	CHRM4 Antibody (BF350)
Species/Host	Rabbit
Reactivity	Human
Conjugation	BF350
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human CHRM4/mAChR M4
Target	CHRM4
Preservatives	1*TBS (pH7.4), 0.05% rAlbumin, 40% Glycerol. Preservative: 0.05% Sodium Azide.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08173
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.