

## Product Datasheet

### CHRM4 Antibody (Cy7) (orb1582148)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | CHRM4 Antibody (Cy7)                                                             |
| <b>Species/Host</b>        | Rabbit                                                                           |
| <b>Reactivity</b>          | Human                                                                            |
| <b>Conjugation</b>         | APC/Cy7                                                                          |
| <b>Tested Applications</b> | FC, ICC                                                                          |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human CHRM4/mAChR M4               |
| <b>Target</b>              | CHRM4                                                                            |
| <b>Preservatives</b>       | 1*TBS (pH7.4), 0.05% rAlbumin, 40% Glycerol. Preservative: 0.05% Sodium Azide.   |
| <b>Form/Appearance</b>     | Liquid                                                                           |
| <b>Concentration</b>       | 1mg/ml                                                                           |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                       |
| <b>Uniprot ID</b>          | <b>P08173</b>                                                                    |
| <b>Dilution Range</b>      | Flow-Cyt=1µg/Test,ICC=1:100-500                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |