


Product Datasheet

PTP zeta antibody (orb158228)

Description	PTP zeta antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human PTP zeta (421-520/2315 aa)
Target	PTP zeta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	251 kDa
Uniprot ID	P23471
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Mouse
spleen...