

Product Datasheet

Pumilio 2 antibody (orb158234)

Description	Pumilio 2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from Human Pumilio 2 (901-1066/1066 aa)
Target	Pumilio 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	114 kDa
Uniprot ID	Q8TB72
Expiration Date	12 months from date of receipt.