

Product Datasheet

rat IgA antibody (orb158264)

Description	rat IgA antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr
Immunogen	Full length plasma protein from rat
Target	rat IgA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	170 kDa
Uniprot ID	P01876
Expiration Date	12 months from date of receipt.