


Product Datasheet

Ron (Phospho-Tyr1238+Tyr1239) antibody (orb158305)

Description	Ron (Phospho-Tyr1238+Tyr1239) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MST1R around the phosphorylation site of Tyr1238+Tyr1239 (E(p-Y)(p-Y)SV)
Target	phospho-Ron (Tyr1238+Tyr1239)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30/119/150 kDa
Uniprot ID	Q04912
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
H-4-II-E
Cel...