



Product Datasheet

Calreticulin Antibody (orb1583181)


Description	Calreticulin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Calreticulin
Target	Calreticulin
Preservatives	1*TBS (pH7.4), 0.05% rAlbumin, 40% Glycerol. Preservative: 0.05% Sodium Azide.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	44 kDa
Uniprot ID	P27797
Dilution Range	WB=1:200-500, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat stomach
using Calr...



IHC-P
analysis of
human brain
using Calr...



IHC-P
analysis of
mouse
stomach
using Ca...