

Product Datasheet

Histone H3 (Di Methyl K4) Antibody (Cy5) (orb1583472)

Description	Histone H3 (Di Methyl K4) Antibody (Cy5)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	APC/Cy5
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of Di Methyl K4
Target	Histone H3 (Di Methyl K4)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	Histone H3.1
Uniprot ID	P68431
Dilution Range	Flow-Cyt=0.2µg /Test,ICC=1:100-500
Expiration Date	12 months from date of receipt.