


Product Datasheet

SNAIL + SLUG(Phospho-Ser246) antibody (orb158432)

Description	SNAIL + SLUG(Phospho-Ser246) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human SNAIL around the phosphorylation site of Ser246 (TF(p-S)RM)
Target	phospho-SNAIL + SLUG (Ser246)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	043623
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Hela Cell
us...