



Product Datasheet

Syntaxin 1A(Phospho-Ser188) antibody (orb158543)

Description	Syntaxin 1A(Phospho-Ser188) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Syntaxin 1a around the phosphorylation site of Ser188 (SI(p-S)KQ)
Target	phospho-Syntaxin 1a (Ser188)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	O95425
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse em...



Immunohistochemical staining of Mouse em...