



Product Datasheet

TIE2 (Phospho-Tyr1108) antibody (orb158606)


Description	TIE2 (Phospho-Tyr1108) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TIE2 around the phosphorylation site of Tyr1108 (TL(p-Y)EK)
Target	Phospho-TIE2 (Tyr1108)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	Q02763
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human
kidney using
Pho...



Western
blot:Lane
1:Mouse
Lung;Lane
2:Mo...



Western
blot:Lane
1:Lung
(Rat);Lane
2:Ki...