

Product Datasheet

**phospho-AKT1 + AKT2 + AKT3
(Thr308+Thr309+Thr305) Antibody (Cy5.5)
(orb1586402)**

Description

phospho-AKT1 + AKT2 + AKT3

(Cy5.5)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

PerCP/Cy5.5

Tested

FC, ICC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human AKT1 around the phosphorylation site of Thr308

Target

phospho-AKT1 + AKT2 + AKT3 (Thr308+Thr309+Thr305)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

P31749

Dilution Range

Flow-Cyt=1µg/Test,ICC=1:100

Expiration Date

12 months from date of receipt.