

## Product Datasheet

**phospho-AKT1 + AKT2 + AKT3  
(Thr308+Thr309+Thr305) Antibody (Cy7)  
(orb1586404)**

**Description**

phospho-AKT1 + AKT2 + AKT3

(Cy7)

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

APC/Cy7

**Tested**

FC, ICC

**Applications**
**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human AKT1 around the phosphorylation site of Thr308

**Target**

phospho-AKT1 + AKT2 + AKT3 (Thr308+Thr309+Thr305)

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Liquid

**Concentration**

1mg/ml

**Storage**

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Uniprot ID**

**P31749**

**Dilution Range**

Flow-Cyt=1µg/Test,ICC=1:100

**Expiration Date**

12 months from date of receipt.