



Product Datasheet

TXNDC11 antibody (orb158674)

Description	TXNDC11 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human TXNDC11 (601-700/985 aa)
Target	TXNDC11
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110 kDa
Uniprot ID	Q6PKC3
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of thyroid ...



Immunohistochemical staining of Rat kidn...