


Product Datasheet

UCMA antibody (orb158683)

Description	UCMA antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human UCMA (21-120/138 aa)
Target	UCMA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14 kDa
Uniprot ID	Q8WVF2
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
young rat
paw using
UC...