

Product Datasheet

**phospho-JNK1 + 2 + 3 (Thr183+Tyr185) Antibody
(BF647) (orb1586854)**

Description

phospho-JNK1 + 2 + 3

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

BF647

Tested

FC, ICC, IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human JNK1 around the phosphorylation site of Thr183/Tyr185

Target

phospho-JNK1 + 2 + 3
(Thr183+Tyr185)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
P45983
Dilution Range

Flow-Cyt=1µg /test,IF=1:100-500,ICC=1:100

Expiration Date

12 months from date of receipt.