


Product Datasheet

UNC80 antibody (orb158689)

Description	UNC80 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human C2orf21 (1901-2100/3258 aa)
Target	UNC80
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	363 kDa
Uniprot ID	Q8N2C7
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...