

Product Datasheet

phospho-ERK1/2 (Thr202 + Tyr204) Antibody (BF350) (orb1586890)

Description

phospho-ERK1/2 (Thr202 + Tyr204)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

BF350

Tested

FC, ICC, IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from mouse p44/42 MAPK around the phosphorylation site of Thr202/Tyr204

Target

phospho-ERK1/2 (Thr202 + Tyr204)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
Q63844
Dilution Range

Flow-Cyt=1µg /test,IF=1:100-500,ICC=1:100

Expiration Date

12 months from date of receipt.