



Product Datasheet

Cadherin 5 antibody (orb158708)


Description	Cadherin 5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human VE Cadherin (251-320/784 aa)
Target	VE Cadherin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81 kDa
Uniprot ID	P33151
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



IHC-P image of human lung tissue using V...



Immunohistochemical staining of human co...



Immunohistochemical staining of paraffin...