

Product Datasheet

VP3 Goose parvovirus antibody (orb158718)

Description	VP3 Goose parvovirus antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from Goose parvovirus VP3 (201-300/534 aa)
Target	Goose parvovirus VP3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	A0A1B2LPK5
Expiration Date	12 months from date of receipt.