


## Product Datasheet

### ZBTB41 antibody (orb158757)

|                            |                                                                                    |
|----------------------------|------------------------------------------------------------------------------------|
| <b>Description</b>         | ZBTB41 antibody                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                       |
| <b>Tested Applications</b> | ELISA, ICC, IF, IHC-Fr, IHC-P, WB                                                  |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human ZBTB41/ZNF924 (701-800/909 aa) |
| <b>Target</b>              | ZBTB41                                                                             |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.              |
| <b>Form/Appearance</b>     | Liquid                                                                             |
| <b>Concentration</b>       | 1mg/ml                                                                             |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.   |
| <b>Note</b>                | For research use only                                                              |
| <b>Isotype</b>             | IgG                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                         |
| <b>MW</b>                  | 105 kDa                                                                            |
| <b>Uniprot ID</b>          | <b>Q5SVQ8</b>                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                    |



Western blot analysis of MCF-7 Human cells...