

Product Datasheet

Phospho-Smad3 (Ser423 + Ser425) Antibody (BF405) (orb1587745)

Description	Phospho-Smad3 (Ser423 + Ser425) Antibody (BF405)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	BF405
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Smad3 around the phosphorylation site of Ser423/425
Target	Phospho-Smad3 (Ser423 + Ser425)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P84022
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.