



Product Datasheet

Laminin alpha 1 antibody (orb1588)


Description	Laminin alpha 1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Laminin alpha 1 (2201-2400/3075 aa)
Target	Laminin alpha 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	337 kDa
Uniprot ID	P25391
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse br...



Immunohistochemical staining of Mouse br...



Immunohistochemical staining of mouse br...