

Product Datasheet

ZNF598(Phospho-Tyr306) antibody (orb158814)

Description	ZNF598(Phospho-Tyr306) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human ZNF598 around the phosphorylation site of Tyr306 (ED(P-Y)EE)
Target	phospho-ZNF598 (Tyr306)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	99 kDa
Uniprot ID	Q86UK7
Expiration Date	12 months from date of receipt.