

## Product Datasheet

### LC3A/B Antibody (BF488) (orb1588233)

|                            |                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------|
| <b>Description</b>         | LC3A/B Antibody (BF488)                                                             |
| <b>Species/Host</b>        | Rabbit                                                                              |
| <b>Reactivity</b>          | Human, Mouse                                                                        |
| <b>Conjugation</b>         | BF488                                                                               |
| <b>Tested Applications</b> | FC, ICC, IF                                                                         |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human LC3A/B                          |
| <b>Target</b>              | LC3A/B                                                                              |
| <b>Preservatives</b>       | Preservative: 15mM Sodium Azide,<br>Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4    |
| <b>Form/Appearance</b>     | Liquid                                                                              |
| <b>Concentration</b>       | 1mg/ml                                                                              |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year.<br>Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                               |
| <b>Isotype</b>             | IgG                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                          |
| <b>Uniprot ID</b>          | <a href="#">Q9GZQ8</a>                                                              |
| <b>Dilution Range</b>      | Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                     |