



Product Datasheet

ZNT8 antibody (orb158829)


Description	ZNT8 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ZNT8. Please contact us for the exact immunogen sequence. The peptide is available as orb374818.
Target	ZNT8
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for one month and at -20°C for a longer time. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	41 kDa
Uniprot ID	Q8BGG0 , P0CE46 , Q8IWU4
Entrez	169026
Dilution Range	IHC-P: 1:100-500, WB: 1:100-500 (based on 0.5mg/ml)
Expiration Date	12 months from date of receipt.



IHC-P image of mouse lymph node tissue u...



Immunohistochemical staining of mouse ly...



IHC-P staining of rat spinal cord tissue...