

Product Datasheet

phospho-c-Jun (Thr239) Antibody (Cy5.5) [Out of stock] (orb1588359)

Description	phospho-c-Jun (Thr239) Antibody (Cy5.5)
Species/Host	Rabbit
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human c-Jun around the phosphorylation site of Thr239
Target	phospho-c-Jun (Thr239)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05412
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.