

Product Datasheet

Phospho-Estrogen Receptor alpha (Ser106) Antibody (Cy5.5) (orb1589988)

Description	Phospho-Estrogen Receptor alpha (Ser106) Antibody (Cy5.5)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	APC/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ER alpha around the phosphorylation site of Ser106
Target	Phospho-Estrogen Receptor alpha (Ser106)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03372
Dilution Range	Flow-Cyt=1µg/Test,ICC=1:100-500
Expiration Date	12 months from date of receipt.