

Product Datasheet

phospho-14-3-3 beta + zeta (Ser186/Ser184) Antibody (PerCP) (orb1591261)

Description	phospho-14-3-3 beta + zeta (Ser186/Ser184) Antibody (PerCP)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PerCP
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human 14-3-3 beta + zeta around the phosphorylation site of Ser186 / Ser184
Target	phospho-14-3-3 beta + zeta (Ser186 / Ser184)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31946
Expiration Date	12 months from date of receipt.