

Product Datasheet

Beta Subunit Cholera Toxin antibody (orb159212)

Description	Beta Subunit Cholera Toxin antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from CTXB/beta subunit Cholera Toxin (51-124/124 aa)
Target	beta subunit Cholera Toxin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	12 kDa
CAS Number	24730-31-2
Uniprot ID	P01556
Expiration Date	12 months from date of receipt.