

Product Datasheet

phospho-Akt (Thr308) Antibody (BF680) (orb1592482)

Description	phospho-Akt (Thr308) Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Akt around the phosphorylation site of Thr308
Target	phospho-Akt (Thr308)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31749
Dilution Range	Flow-Cyt=1µg /test,IF=1:100-500
Expiration Date	12 months from date of receipt.