



Product Datasheet

Cyp2S1 antibody (orb159252)

Description	Cyp2S1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CYP2S1 (21-120/504 aa)
Target	CYP2S1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	Q96SQ9
Dilution Range	IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, WB=1:500-2000, ICC=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human liver tissue.



Immunohistochemical staining of Human gastric tissue.