

## Product Datasheet

### Insulin Receptor Beta Antibody (BF680) (orb1592580)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | Insulin Receptor Beta Antibody (BF680)                                           |
| <b>Species/Host</b>        | Rabbit                                                                           |
| <b>Reactivity</b>          | Human                                                                            |
| <b>Conjugation</b>         | BF680                                                                            |
| <b>Tested Applications</b> | FC, ICC, IF                                                                      |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human Insulin Receptor Beta        |
| <b>Target</b>              | Insulin Receptor Beta                                                            |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.            |
| <b>Form/Appearance</b>     | Liquid                                                                           |
| <b>Concentration</b>       | 1mg/ml                                                                           |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                       |
| <b>CAS Number</b>          | 3                                                                                |
| <b>Uniprot ID</b>          | <b>P06213</b>                                                                    |
| <b>Dilution Range</b>      | Flow-Cyt=0.2µg/Test,IF=1:100-500,ICC=1:100                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |