


Product Datasheet

IFT46 antibody (orb159283)

Description	IFT46 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human IFT46 (101-200/304 aa)
Target	IFT46
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34 kDa
Uniprot ID	Q9NQC8
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Molt-4(Human) Cell; K562(Hum...)