

Product Datasheet

KLF4 Antibody (BF750) (orb1593238)

Description	KLF4 Antibody (BF750)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	BF750
Tested Applications	FC, ICC, IF
Immunogen	Recombinant protein within human KLF4 aa 1-500
Target	KLF4
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	O43474
Dilution Range	Flow-Cyt=1ug/Test,IF=1:50-200,ICC=1:100
Expiration Date	12 months from date of receipt.