

Product Datasheet

KLF4 Antibody (HRP) (orb1593255)

Description	KLF4 Antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	Recombinant protein within human KLF4 aa 1-500
Target	KLF4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	O43474
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:50-200,ICC=1:100
Expiration Date	12 months from date of receipt.